















VISIBLE METAL DESIGN CHART

FOR PORCELAIN TO METAL RESTORATIONS

This chart has been devised as a time-saving aid in prescribing conventional porcelain to metal restorations. Using these drawings as your guide should also reduce or eliminate the need for remakes. If you find that slight modifications in this chart will suit your individual preference, please let us know. We want to make working together as easy and efficient as possible.

INSTRUCTIONS:

1. Both the dentist and technician have identical charts.
2. The dentist clarifies Metal Design for each case by simply writing the corresponding letter on the prescription, avoiding time-consuming long hand instructions.
3. The technician constructs the case according to the dentist's designated letter, eliminating misinterpretation.

 <p>A</p> <p>Full Porcelain, no metal showing</p>	 <p>B</p> <p>Lingual Shoulder</p>	 <p>C</p> <p>Lingual shoulder with small labial shoulder</p>	 <p>D</p> <p>Narrow shoulder all around</p>	 <p>E</p> <p>Full metal lingual</p>	 <p>F</p> <p>Full metal lingual with narrow labial shoulder</p>
 <p>G</p> <p>Full Porcelain, no metal showing</p>	 <p>H</p> <p>Lingual Shoulder</p>	 <p>J</p> <p>Lingual shoulder with small labial shoulder</p>	 <p>K</p> <p>Lingual cusps metal</p>		
 <p>L</p> <p>Lingual cusps metal with narrow buccal shoulder</p>	 <p>M</p> <p>Buccal porcelain veneer</p>	 <p>N</p> <p>Buccal porcelain veneer with narrow buccal shoulder</p>	 <p>P</p> <p>Buccal porcelain veneer with protected buccal cusp</p>		